




University
of Southampton



THE ELEMENTS OF SUCCESS

Ionisation in the absence of high voltage using SFC-MS: a route to greater sensitivity

Mohini Thite

Email: mat@soton.ac.uk


Outline

- What is SFC?
- Coupling SFC-MS
- Experimental Design
- Detection
 - ESI-APCI
 - SFC conditions
 - Pressure
 - Modifier
 - Flow rate

What is SFC?

- Supercritical fluid chromatography
 - Mobile phase in supercritical state
- Mobile phase: non polar carbon dioxide
- Modifier: polar organic solvent to improve chromatography
- Polar stationary phases
 - e.g. silica, amine, 2-EP
- Considered as normal phase technique

SFC-MS Coupling



Ref: M. Garzotti, M. Hamdan, *J. Chromatogr B*, 770 (2002) 53-61.

MS Detection

- Initial studies utilised ESI
- Optimisation of API for SFC
- Observation of samples ionising in the absence of high voltages
- Specific test compounds analysed to probe this ionisation mechanism


What is causing ionisation?

- API source configuration dependent?
 - ESI *cf.* APCI
- Chromatographic system specific?
 - HPLC *cf.* SFC

Is ionisation related to?

- Charge residual model?

..... r . i




• /% molecule charged

- Nebuliser gas


Ref: C. Dass, Principles and Practice of Biological Mass Spectrometry, John Wiley & Sons, Inc., New York, 2001.

Experimental conditions





	1	2	3	4
Pure methanol	HV: 3.5kV CV: 20V	HV: 3.5kV CV: 0V	HV: 0kV CV: 20V	HV: 0kV CV: 0V
Methanol + 0.1% NH ₄ OAc	HV: 3.5kV CV: 20V	HV: 3.5kV CV: 0V	HV: 0kV CV: 20V	HV: 0kV CV: 0V
Methanol + 0.1% HCOOH	HV: 3.5kV CV: 20V	HV: 3.5kV CV: 0V	HV: 0kV CV: 20V	HV: 0kV CV: 0V
No make-up	HV: 3.5kV CV: 20V	HV: 3.5kV CV: 0V	HV: 0kV CV: 20V	HV: 0kV CV: 0V

ESI reserpine HV on/off






ESI and APCI



Summary of ESI and APCI


- Comparable/Improved ionisation in absence of high voltage
- Ionisation suppressed with addition of ammonium acetate
- Formic acid does not enhance ionisation
- Data for APCI source comparable with ESI source
- Ionisation independent of capillary type

Pressure and modifier




- Mixture of terfenadine, oxybutynin and erythromycin analysed at varying pressure and modifier percentage
- Optimum conditions for maximum sensitivity in presence or absence of high voltage:
Pressure 100bar and modifier 30% MeOH

Flow rate




High flow rate (5mL/min)




- Improvement in sensitivity upon increasing the distance between the ESI needle tip and the sample cone

Modifier effect

Comparison of acetonitrile and methanol modifiers




Modifier effect



- Increased ionisation with MeOH modifier compared to MeCN
- Is the ionisation efficiency related to preformed ions in solution?

High voltage profile



- Not electrospray
- CRM - Ions formed in solution or during spray process
- CRM plus other processes, sonic spray


Conclusions

- Compounds ionise in the absence of high voltage
∴ not electrospray ionisation
- Addition of ammonium acetate does not aid
ionisation ∴ not thermospray ionisation
- Suggests ions are formed in solution or during the
supersonic jet expansion

Conclusions

- MeCN as make-up solvent or as modifier does not assist in ionisation, unlike MeOH ∴ supporting CRM?
- Ionisation on FIA-ESI-MS with HV turned off but not as good as SFC-MS
- API-MS response dependent on modifier concentration; this should be taken into account for gradient elution studies

Conclusions



For small molecule analysis if sensitivity is an issue :
Turning off the electrospray high voltage could
improve LOD - e.g. oxybutynin, reserpine data

Acknowledgements

G. John Langley
Julie Herniman
Amaury Cazenave Gassiot

School of Chemistry,
University of Southampton,
SO17 1BJ, UK

Paul Oakley

Mettler-Toledo Autochem, 130
Executive Drive, Suite 2A, Newark,
DE 19702, USA

Laure Hitzel
Frank Pullen

Pfizer Global Research and
Development, Sandwich,
CT13 9NJ, UK

Jeff Caldwell

Princeton Chromatography Inc.,
Cranbury, NJ 08512, USA

Clare Paterson
Bob Boughtflower

GlaxoSmithKline, CASS, New
Frontiers North, Harlow,
CM19 5AW, UK

Steve Lane

GlaxoSmithKline, Gunnels Wood
Road, Stevenage, SG1 2NY, UK

Stefan Richardson

Evotec, 151 Milton
Park, Abingdon,
OX14 4SD, UK



BERGER SFC™



University
of Southampton



BRITISH MASS SPECTROMETRY SOCIETY



evotec



GlaxoSmithKline

Thank You